DIGITHEAT 19L	2001241
DIGITHEAT 36L	2001242
DIGITHEAT 52L	2001243
DIGITHEAT 80L	2001244
DIGITHEAT 150L	2001245

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***DRYING AND STERILIZING OVENS NATURAL CONVECTION
WITH TEMPERATURE AND TIME REGULATION
AND DIGITAL READING***

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GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) *Do not use the apparatus with liquids which can give off vapours capable of making explosive mixtures.*

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CONTENTS

GENERAL INFORMATION
CONTENTS
PACKING LISTS
ACCESSORIES
EQUIPMENT DESCRIPTION
TECHNICAL FEATURES
INSTALLATION
OPERATION
MAINTENANCE
SPARE PARTS
ELECTRICAL DIAGRAM
GUARANTEE
"EC" CONFORMITY DECLARATION

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The standard equipment consist of the following components:

Pieza / Piece	Código bandeja Shelf code	Cant. / Qty.	Código guías Shelf runners code	Cant. / Qty.
Estufa / Oven 2001241	2000021	2	2000011	4
Estufa / Oven 2001242	2000022	2	2000012	4
Estufa / Oven 2001243	2000024	2	2000014	4
Estufa / Oven 2001244	2000023	2	2000013	4
Estufa / Oven 2001245	2000025	2	2000015	4
Manual de instrucciones Instruction Manual	80034	1	--	--

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ACCESSORIES

(The following components must be installed in factory)

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SAFETY THERMOSTAT <i>Installed in the control panel, adjustable, it disconnects the heating in case of oven's own thermostat failure.</i>	2000001
--	---------

TEMPERATURE CHART RECORDER <i>With graphic print on continuous paper.</i>	2000019
--	---------

DIGITAL TEMPERATURE RECORDER <i>With numerical print on continuous paper with print intervals from 1" to 99 hours.</i>	2000016
---	---------

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EQUIPMENT DESCRIPTION

DIGITHEAT ovens are designed for drying and sterilization for adjustable temperatures from ambient +5 °C up to 250°C.

The oven is provided with a digital electronic temperature and time regulator, and with a Pt100 probe.

The heat distribution inside the chamber is produced by natural convection.

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Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

	2001241	2001242	2001243	2001244	2001244	
Capacidad / capacity (L.)	19	36	52	80	150	
Temperatura máxima / Maximum temperature	250°C					
*Homogeneidad / Homogeneity	± 2%					
*Estabilidad / Stability	± 0,25 °C					
Error de consigna / Thermometer error	± 2%					
Resolución / Resolution	1°C					
Medidas interiores Internal dimensions cm.	Alto / Height	30	40	33	50	50
	Ancho / Length	25	30	47	40	60
	Fondo / Width	25	30	33	40	50
Medidas exteriores Overall dimensions cm.	Alto / Height	50	60	53	70	70
	Ancho / Length	60	65	82	74	95
	Fondo / Width	44	49	52	59	68
Consumo / Consumption (W.)	600	900	1000	1200	2000	
Peso / Weight (Kg.)	24	33	40	51	68	

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** Measured according to the measurement of J.P. Selecta, s.a.*

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INSTALLATION

Put the oven on a flat, horizontal, level surface leaving a space of about 10 cm at the back of the apparatus. For a better oven function, do not situate it in places where it can be directly exposed to the sun radiation or heat sources.

For an optimum function the ambient temperature must not be over 30°C.

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Do not use the oven to dry or for thermal treatment of products which can give off vapours capable of making an explosive mixture.

Be sure that the voltage supply is the same as the one indicated on the characteristics plate.

Do not use the apparatus if it is not earthed.

If you change the plug bear the following in mind:


**Blue cable: Neutral.
Brown cable: Phase.
Yellow/green cable: Earth.**

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CONTROL PANEL:

2. *Temperature mode indicator.*
3. *Time mode indicator.*
4. *Temperature/Time display.*
5. *Starting up indicator lamp.*
6. *Waiting time indicator lamp.*
7. *Temperature selector push-button (Set of temperature from ambient +5°C up to 250°C).*
8. *Time selector push-button.*
 - *Set of starting up time: Working period between 1' and 9h 59', or up to 99.9h. once the set temperature has been reached.*
 - *Set of waiting time: Period of time until the starting up, between 1 and 24 hours. Decimal fraction of hour (1.5 = 1hour 30 minutes).*
9. *Parameter increase push-button.*
10. *Parameter decrease push-button.*
11. *Start/Stop push-button.*

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The oven can work in 4 ways:

1. NORMAL MODE:

When pressing the main switch (1), placed in the lower part of the frontal, the display lights up (4).

Pressing the push-button (7) and the (9) or (10) we preset the desired temperature.

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To stop, press the push-button again (11) and the indicator will be turned off (5).

2. TEMPORIZED MODE AT THE DISCONNECTION:

When pressing the main switch (1) the display lights up (4).

Pressing simultaneously the push-button (7) and the (9) or (10) we preset the desired temperature.

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Pressing simultaneously the push-button (8) and the (9) or (10) we preset the disconnection time.

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The disconnection time begins when the oven has arrived to the preset temperature.

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In this operation mode, the indicator (5) remains intermittent until the oven reaches to the preset temperature, once reached, the indicator (5) remains lighted up.

When the preset time has lapsed, the oven stops automatically and the indicator (5) will go off.

You can stop it manually pressing the push-button (11) again.

3. TEMPORIZED MODE AT THE STARTING UP:

When pressing the main switch (1) the display lights up (4).

Pressing simultaneously the push-button (7) and the (9) or (10) we preset the desired temperature.

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Pressing simultaneously the push-button (8) twice and the (9) or (10) we preset the waiting time until the starting up. During this operation the indicator (6) will light up.

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In this operation mode, the indicator (6) remains lit up during the waiting time, this over, the indicator lamp (6) goes off and the oven starts automatically.

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4. TEMPORIZED WAY AT THE STARTING UP AND AT THE DISCONNECTION:

When pressing the main switch (1) the display lights up (4).

Simultaneously pressing the push-button (7) and the (9) or (10) we preset the desired temperature.

During this operation the indicator (2) lights up.

Simultaneously pressing the push-button (8) twice and the (9) or (10) we preset the waiting time until the starting up. During this operation the indicator lamp (6) will light up.

Simultaneously pressing the push-button (8) and the (9) or (10) we preset the disconnection time.

During this operation the indicator (3) lights up.

The disconnection time begins when the oven has reached to the preset temperature.

Pressing the push-button (11) the oven will start and the indicator (5) will light up.

In this operation mode, the indicator lamp (5) remains intermittent until the oven reaches to the preset temperature, once reached, the indicator (5) remains lighted up.

Lapsed the preset time, the oven stops automatically and the indicator (5) will go off.

You can stop it manually by pressing the push-button (11) again.

Pressing the push-button (11) the oven will start and the indicators (5) and (6) will light up.

In this operation mode, the indicator lamp (6) remains lit up during the waiting time, once over, the indicator lamp (6) goes off and the oven starts automatically.

OVENS WITH SAFETY THERMOSTAT:

To adjust the safety temperature thermostat do as follows:

1. *With the oven connected to the mains, turn the safety thermostat control knob clockwise to the end of the scale.*
2. *Ensure that the safety thermostat is reset by pressing the reset button. The red indicator lamp must be off.*
3. *Select the desired safety temperature with the temperature control knob. This temperature has to be minimum at 5°C above the working temperature.*
4. *With the door closed, turn on the oven with the switch.*

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5. *When the temperature of the oven is stabilized, turn the safety thermostat control knob slowly to the left until it clicks, the red indicator lamp lights up and the heating turns off permanently.*
6. *Stop the oven with the switch. Open the door and leave it to cool with the door open until the temperature goes down by 50% of the working temperature (approx.).*
7. *Reset the safety thermostat by pressing the reset button, the red indicator lamp must turn off.*
8. *Now the oven is ready to work.*

NOTE: *The safety thermostat is useful when the user works habitually at the same temperature.*

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MAINTENANCE

Before removing the casing disconnect the apparatus from the mains.

The manipulation of the internal electronic circuits of the oven by unauthorized personnel can cause irreparable damage. Take it to one of the J.P.SELECTA authorized technical services.

CLEANING:

For the cleaning of the different parts use the following products:

Cleaning of stainless steel: Alcohol.

Cleaning of plastic: Alcohol with cotton duster.

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
SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P.SELECTA,S.A.

Descripción / Description	Código / Code				
	Modelo / Model	2001241	2001242	2001243	2001244
Junta silicona / Silicone gasket	21067	21068	21077	21069	21078
Solid relé / Relay	36026	36027			36022
Sonda Pt100 / Pt100 probe	43074				
Resistencia / Heating element	39122	39118	39120	39119	39121
Cierre / Lock	10012				
Circuito electrónico / P.C.B	29130				

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GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA, S.A.

Any manipulation of the apparatus by personnel not authorized by J.P. SELECTA, S.A. automatically cancels the guarantee.

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